## 250 $\mu\text{L},$ Model 1725 AD SYR, M8 Threads

250 µL, Model 1725 AD SYR, M8 Threads



Part Number/REF: 201100 Unit Price: Quote Required

## Specifications

Needle	No Needle Available
Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	AccuDil® (AD)
Std. Needle Dead Volume	N/A
Autoclavable	No
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.