$500~\mu\text{L},\,\text{Model}\,\,1750~\text{AD}\,\,\text{SYR},\,\text{M8}\,\,\text{Threads}$

 $500~\mu\text{L},$ Model 1750 AD SYR, M8 Threads



Part Number/REF: 201150 Unit Price: Quote Required

Specifications

Needle	No Needle Available
Volume	500 μL
Syringe Type	Gastight
Syringe Series	1700
Termination	AccuDil® (AD)
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	700 psig (48.3 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.

1 of 1 1/4/2017 2:20 PM