

# 1 L, Model S1000 TLL SYR, NDL Sold Separately

1 L, Model S1000 TLL SYR, NDL Sold Separately



Part Number/REF: 86312

Unit Price: Quote Required

## Specifications

|                         |   |
|-------------------------|---|
| Volume                  | 1 L   |
| Syringe Type            | Gastight                                    |
| Syringe Series          | Super Syringe                               |
| Termination             | PTFE Luer Lock (TLL)                        |
| Needle                  | Metal (N) or Kel-F                          |
| Gauge                   | No Needle Included                          |
| Point Style             | No Needle Included                          |
| Needle Length           | No Needle Included                          |
| Std. Needle Dead Volume | N/A   |
| Fluid Path              | Buna-N, Nickel plated brass, PTFE, Silicone |
| Autoclavable            | No  |
| Gas Sterilizable        | Yes   |
| Pressure                | 30 psig (2.1 bar)                           |
| Maximum Temperature     | 80°C / 176°F                                |
| Minimum Temperature     | 0°C / 32°F                                  |
| Barrel Inner Diameter   | 3.44 inches (87.4 mm)                       |
| Barrel Outer Diameter   | 3.95 inches (100 mm)                        |
| Sales Quantity          | 1 syringe per pack                          |

## Accessories

| Part Number/REF | Description                           | Unit Price     |
|-----------------|---------------------------------------|----------------|
| 18026           | TLL Connector S0500 S1000 S1500 S2000 | Quote Required |

