

# 5 mL, Model 1005 AD SYR, M8 Threads

---

5 mL, Model 1005 AD SYR, M8 Threads

Part Number/REF: 201300

Unit Price: Quote Required

## Specifications

---

Needle	No Needle Available
Volume	5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	AccuDil® (AD)
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

---

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.