

500 μ L, Model 1750 Terasaki Dispenser

500 μ L, Model 1750 Terasaki Dispenser

Part Number/REF: 201940

Unit Price: Quote Required

Specifications

Volume	500 μ L
Syringe Type	Gastight
Syringe Series	1700
Fluid Path	Borosilicate glass, PTFE
Sales Quantity	1 dispenser

Accessories

Part Number/REF	Description	Unit Price
10370A	Glass Barrel 1750	Quote Required
1169-01	Plunger Assembly N RN 1750	Quote Required

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.