

25 mL, Model 1025 AD SYR, M8 Threads

25 mL, Model 1025 AD SYR, M8 Threads

Part Number/REF: 201400

Unit Price: Quote Required

Specifications

Needle	No Needle Available
Volume	25 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	AccuDil® (AD)
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
13271	Plunger Assembly LTLNTLL w wo Slots ADDAD 1025	Quote Required

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.