

500 μ L, Model 1750 DAD SYR, Diluter Syringe

500 μ L, Model 1750 DAD SYR, Diluter Syringe

Part Number/REF: 201650

Unit Price: Quote Required

Specifications

Autoclavable	No
Needle	No Needle Available
Volume	500 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	Diluter AccuDil® (DAD)
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	20°C / 68°F
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.