

# 25 $\mu$ L, Model 1702 Terasaki Dispenser

---

25  $\mu$ L, Model 1702 Terasaki Dispenser

Part Number/REF: 201900

Unit Price: Quote Required

## Specifications

---

Volume	25 $\mu$ L
Syringe Type	Gastight
Syringe Series	1700
Fluid Path	Borosilicate glass, PTFE
Sales Quantity	1 dispenser

---

## Accessories

---

Part Number/REF	Description	Unit Price
1122-01	1702 Plunger Assembly LT N RN C	Quote Required

---

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.