

500 µL, Model 1750 AD-SAL SYR, M8 Threads

500 µL, Model 1750 AD-SAL SYR, M8 Threads

Part Number/REF: 203320

Unit Price: Quote Required

Specifications

Needle	No Needle Available
Volume	500 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	AccuDil® (AD)
Fluid Path	Borosilicate glass, PTFE
Plunger Style	Saltline (SAL)
Autoclavable	No
Gas Sterilizable	Yes
Pressure	700 psig (48.3 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Outer Diameter	0.028 inch (0.7 mm)
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.