

# 1 mL, Model 1001 C SYR, 1/4-28 Threads

---

1 mL, Model 1001 C SYR, 1/4-28 Threads

Part Number/REF: 203540

Unit Price: Quote Required

## Specifications

Volume	1 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	ChemSeal (C)
Needle	No Needle Available
Fluid Path	Borosilicate glass, PTFE
Plunger Style	SaltLine w/ Cavo 30 mm Plunger Button (XP)
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Outer Diameter	0.028 inch (0.7 mm)
Sales Quantity	1 syringe per pack

---

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.