

# 5 mL, Model 1005 C-XP, Instrument Syringe

---

5 mL, Model 1005 C-XP, Instrument Syringe

Part Number/REF: 208645

Unit Price: Quote Required

## Specifications

Maximum Temperature	122°F (50°C)
Minimum Temperature	68°F (20°C)
Termination	ChemSeal (C)
Fluid Path	Borosilicate glass, PTFE
Sales Quantity	1 syringe per pack
Volume	5 mL
Gas Sterilizable	Yes
Syringe Series	1000
Maximum Pressure	200 psig (13.8 bar)
Autoclavable	Yes
Syringe Type	Gastight
Needle	No Needle Available

---

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.