

2.5 mL, Model 1002 RN SYR, NDL Sold Separately

2.5 mL, Model 1002 RN SYR, NDL Sold Separately



Part Number/REF: 7650-01

Unit Price: Quote Required

Specifications

Volume	2.5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.287 inches (7.29 mm)
Barrel Outer Diameter	0.407 inches (10.3 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
1360-01	Plunger Assembly LTLTNTLL w amp wo Slots C 1002	Quote Required
30902	RN Nut	Quote Required
CAL7650-01	25 mL Model 1002 RN Calibrated SYR NDL Sold Separately	Quote Required

Accessories

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.