

5 mL, Model 1005 RN SYR, NDL Sold Separately

5 mL, Model 1005 RN SYR, NDL Sold Separately



Part Number/REF: 7651-01

Unit Price: Quote Required

Specifications

| | |
|-------------------------|--------------------------|
| Volume | 5 mL |
| Syringe Type | Gastight |
| Syringe Series | 1000 |
| Termination | Removable Needle (RN) |
| Needle | Large RN |
| Gauge | No Needle Included |
| Point Style | No Needle Included |
| Needle Length | No Needle Included |
| Std. Needle Dead Volume | N/A |
| Fluid Path | Borosilicate glass, PTFE |
| Autoclavable | Yes |
| Gas Sterilizable | Yes |
| Pressure | 200 psig (13.8 bar) |
| Maximum Temperature | 115°C / 239°F |
| Minimum Temperature | 10°C / 50°F |
| Barrel Inner Diameter | 0.406 inches (10.3 mm) |
| Barrel Outer Diameter | 0.530 inches (13.5 mm) |
| Sales Quantity | 1 syringe per pack |

Accessories

| Part Number/REF | Description | Unit Price |
|-----------------|---|----------------|
| 13230 | Plunger Assembly LTLTNTLL w wo Slots RN 1005 | Quote Required |
| 30902 | RN Nut | Quote Required |
| CAL7651-01 | 5 mL Model 1005 RN Calibrated SYR NDL Sold Separately | Quote Required |

Accessories

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.