

# 10 mL, Model 1010 RN SYR, NDL Sold Separately

10 mL, Model 1010 RN SYR, NDL Sold Separately



Part Number/REF: 7652-01

Unit Price: Quote Required

## Specifications

Volume	10 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.574 inches (14.6 mm)
Barrel Outer Diameter	0.695 inches (17.7 mm)
Sales Quantity	1 syringe per pack

## Accessories

Part Number/REF	Description	Unit Price
13231	Plunger Assembly LTLTNTLL w wo Slots RN 1010	Quote Required
30902	RN Nut	Quote Required
CAL7652-01	10 mL Model 1010 RN Calibrated SYR NDL Sold Separately	Quote Required

## Accessories

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.