


10 µL, Model 1701 RN SYR, NDL Sold Separately

10 µL, Model 1701 RN SYR, NDL Sold Separately



Part Number/REF: 7653-01

Unit Price: Quote Required

Specifications

Volume	10 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.018 inches (0.460 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
13205	Plunger Assembly N, RN, LT, LTN 1701	Quote Required
30902	RN Nut	Quote Required
CAL7653-01	10 microL Model 1701 RN Calibrated SYR NDL Sold Separately	Quote Required

Accessories

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.