

# 100 $\mu$ L, Model 1710 RN SYR, NDL Sold Separately

100  $\mu$ L, Model 1710 RN SYR, NDL Sold Separately



Part Number/REF: 7656-01

Unit Price: Quote Required

## Specifications

Volume	100 $\mu$ L
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

## Accessories

Part Number/REF	Description	Unit Price
1162-02	Plunger Assembly N RN LTLN SampleLock TLL w wo Slots 1710	Quote Required
30902	RN Nut	Quote Required
CAL7656-01	100 microL Model 1710 RN Calibrated SYR NDL Sold Separately	Quote Required

## Accessories

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.