

10 μ L, Model 1701 LT SYR, NDL Sold Separately

10 μ L, Model 1701 LT SYR, NDL Sold Separately



Part Number/REF: 80001

Unit Price: Quote Required

Specifications

Volume	10 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	Luer Tip (LT)
Needle	Kel-F Hub
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Autoclavable	Yes to 115°C
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.018 inches (0.460 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
13205	Plunger Assembly N, RN, LT, LTN 1701	Quote Required
CAL80001	10 microL Model 1701 LT Calibrated SYR NDL Sold Separately	Quote Required

Accessories

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.