

10 µL, Model 1701 RN SYR, Small Removable NDL, 32 ga, 2 in, point style 3

10 µL, Model 1701 RN SYR, Small Removable NDL, 32 ga, 2 in, point style 3



Part Number/REF: 80014
Unit Price: Quote Required

Specifications

Volume	10 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	32 gauge
Point Style	3
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0045 inches (0.11 mm)
Needle Outer Diameter	0.00925 inches (0.24 mm)
Std. Needle Dead Volume	0.47 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.018 inches (0.460 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7762-05	32 gauge Small Hub RN NDL 2 in point style 3 6PK	Quote Required

13205	Plunger Assembly N, RN, LT, LTN 1701	Quote Required
30902	RN Nut	Quote Required

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.