

10 μ L, Model 1701 CX SYR, 1/4-28 Threads

10 μ L, Model 1701 CX SYR, 1/4-28 Threads



Part Number/REF: 80062

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	10 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	ChemSeal (C)
Needle	No Needle Available
Plunger Style	Instrument Plunger Stop (X)
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.018 inches (0.460 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.