

50 μ L, Model 1705 SL SYR, Large Removable NDL, 22s ga, 2 in, point style 2

50 μ L, Model 1705 SL SYR, Large Removable NDL, 22s ga, 2 in, point style 2



Part Number/REF: 80956

Unit Price: Quote Required

Specifications

Volume	50 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	Sample Lock (SL)
Needle	Large RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0059 inches (0.15 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	1.13 μ L
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.041 inches (1.03 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7779-03	22s gauge Large Hub RN NDL 2 in point style 2 6PK	Quote Required
1162-01	Plunger Assembly N RN LTLN SampleLock TLL w wo Slots 1705	Quote Required
30902	RN Nut	Quote Required

Accessories



Large Hub Removable Needles

These needles are available in 20-26 gauges. The standard length is 2 inches. Custom lengths are from 0.375 to 12 inches.

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.