

50 µL, Model 1705 CX SYR, 1/4-28 Threads

50 µL, Model 1705 CX SYR, 1/4-28 Threads



Part Number/REF: 80962

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	50 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	ChemSeal (C)
Needle	No Needle Available
Plunger Style	Instrument Plunger Stop (X)
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.041 inches (1.03 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
1117-01	Plunger Assembly TLLXDXCX 1705	Quote Required