

100 μ L, Model 1710 CX SYR, 1/4-28 Threads

100 μ L, Model 1710 CX SYR, 1/4-28 Threads



Part Number/REF: 81060

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	100 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	ChemSeal (C)
Needle	No Needle Available
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.