

100 μ L, Model 1710 CX SYR, 1/4-28 Threads, Plunger

100 μ L, Model 1710 CX SYR, 1/4-28 Threads, Plunger



Part Number/REF: 81062

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	100 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	ChemSeal (C)
Needle	No Needle Available
Plunger Style	Instrument Plunger Stop (X)
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
1117-02	Plunger Assembly TLLXDXCX 1710	Quote Required