

250 µL, Model 1725 DX SYR, Diluter Syringe

250 µL, Model 1725 DX SYR, Diluter Syringe



Part Number/REF: 81126

Unit Price: Quote Required

Specifications

Autoclavable	No
Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Diluter PTFE Luer Lock (D)
Needle	No Needle Available
Plunger Style	Instrument Plunger Stop (X)
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
1117-03	Plunger Assembly TLLDXCX 1725	Quote Required