

250 µL, Model 1725 SL SYR, Large Removable NDL, 22s ga, 2 in, point style 2

250 µL, Model 1725 SL SYR, Large Removable NDL, 22s ga, 2 in, point style 2



Part Number/REF: 81156
Unit Price: Quote Required

Specifications

Volume	250 µL
Syringe Type	Gastight
Syringe Series	1700
Termination	Sample Lock (SL)
Needle	Large RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0059 inches (0.15 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	1.13 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7779-03	22s gauge Large Hub RN NDL 2 in point style 2 6PK	Quote Required

1320-01	Plunger Assembly SampleLock 1725	Quote Required
30902	RN Nut	Quote Required
CAL81156	250 microL Model 1725 SL Calibrated SYR Large Removable NDL 22s ga 2 in point style 2	Quote Required

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.