

250 μ L, Model 1725 C SYR, 1/4-28 Threads

250 μ L, Model 1725 C SYR, 1/4-28 Threads



Part Number/REF: 81160

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	250 μ L
Syringe Type	Gastight
Syringe Series	1700
Termination	ChemSeal (C)
Needle	No Needle Available
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.