

# 500 $\mu$ L, Model 1750 CX SYR, 1/4-28 Threads

500  $\mu$ L, Model 1750 CX SYR, 1/4-28 Threads



Part Number/REF: 81262

Unit Price: Quote Required

## Specifications

|                       |                               |
|-----------------------|-------------------------------|
| Fluid Path            | Borosilicate glass, PTFE, PFA |
| Volume                | 500 $\mu$ L                   |
| Syringe Type          | Gastight                      |
| Syringe Series        | 1700                          |
| Termination           | ChemSeal (C)                  |
| Needle                | No Needle Available           |
| Plunger Style         | Instrument Plunger Stop (X)   |
| Autoclavable          | Yes                           |
| Gas Sterilizable      | Yes                           |
| Pressure              | 700 psig (48.3 bar)           |
| Maximum Temperature   | 115°C / 239°F                 |
| Minimum Temperature   | 10°C / 50°F                   |
| Barrel Inner Diameter | 0.128 inches (3.26 mm)        |
| Barrel Outer Diameter | 0.305 inches (7.75 mm)        |
| Sales Quantity        | 1 syringe per pack            |

## Accessories

| Part Number/REF | Description                    | Unit Price     |
|-----------------|--------------------------------|----------------|
| 1120-01         | Plunger Assembly TLLXDXCX 1750 | Quote Required |