

1 mL, Model 1001 LT SYR, NDL Sold Separately

1 mL, Model 1001 LT SYR, NDL Sold Separately



Part Number/REF: 81301

Unit Price: Quote Required

Specifications

Volume	1 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip (LT)
Needle	Kel-F Hub
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Autoclavable	Yes to 115°C
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.181 inches (4.61 mm)
Barrel Outer Diameter	0.355 inches (9.02 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
1359-01	Plunger Assembly LTLTNTLL w wo Slots SampleLock C 1001	Quote Required
CAL81301	1 mL Model 1001 LT Calibrated SYR NDL Sold Separately	Quote Required

Accessories

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.