

1 mL, Model 1001 DX SYR, Diluter Syringe

1 mL, Model 1001 DX SYR, Diluter Syringe



Part Number/REF: 81326

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	1 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Diluter PTFE Luer Lock (D)
Needle	No Needle Available
Plunger Style	Instrument Plunger Stop (X)
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.181 inches (4.61 mm)
Barrel Outer Diameter	0.355 inches (9.02 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
1578-01	Plunger Assembly DX 1001	Quote Required
CAL81326	1 mL Model 1001 DX Calibrated SYR Diluter Syringe	Quote Required