

1 mL, Model 1001 RN SYR, Large Removable NDL, 22 ga, 2 in, point style 2

1 mL, Model 1001 RN SYR, Large Removable NDL, 22 ga, 2 in, point style 2



Part Number/REF: 81330
Unit Price: Quote Required

Specifications

Volume	1 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	22 gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	6.81 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.181 inches (4.61 mm)
Barrel Outer Diameter	0.355 inches (9.02 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7779-01	22 gauge Large Hub RN NDL 2 in point style 2 6PK	Quote Required
1359-01	Plunger Assembly LTLTNTLL w wo Slots SampleLock C 1001	Quote Required
30902	RN Nut	Quote Required

CAL81330

1 mL Model 1001 RN Calibrated SYR Large Removable NDL 22 ga 2 in point style 2

Quote Required

Accessories



PB600-1 Repeating Dispenser

The PB600-1 is a sampling device designed for repeated dispensing of equal volumes. The instrument delivers 1/50th of the syringe volume.



Large Hub Removable Needles

These needles are available in 20-26 gauges. The standard length is 2 inches. Custom lengths are from 0.375 to 12 inches.

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.