

# 2.5 mL, Model 1002 LTN SYR, Cemented NDL, 22 ga, 2 in, point style 3

2.5 mL, Model 1002 LTN SYR, Cemented NDL, 22 ga, 2 in, point style 3



Part Number/REF: 81416

Unit Price: Quote Required

## Specifications

Volume	2.5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip Cemented Needle (LTN)
Needle	Cemented Needle
Gauge	22 gauge
Point Style	3
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	6.81 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.287 inches (7.29 mm)
Barrel Outer Diameter	0.407 inches (10.3 mm)
Sales Quantity	1 syringe per pack

## Accessories

Part Number/REF	Description	Unit Price
1360-01	Plunger Assembly LTLTNTLL w amp wo Slots C 1002	Quote Required
CAL81416	25 mL Model 1002 LTN Calibrated SYR Cemented NDL 22 ga 2 in point style 3	Quote Required

## Accessories



### PB600-1 Repeating Dispenser

The PB600-1 is a sampling device designed for repeated dispensing of equal volumes. The instrument delivers 1/50th of the syringe volume.

---

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.