

5 mL, Model 1005 TLL SYR, NDL Sold Separately

5 mL, Model 1005 TLL SYR, NDL Sold Separately



Part Number/REF: 81520

Unit Price: Quote Required

Specifications

Autoclavable	Yes
Volume	5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Needle	Metal (N) or Kel-F
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.406 inches (10.3 mm)
Barrel Outer Diameter	0.530 inches (13.5 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
13230	Plunger Assembly LTLTNTLL w wo Slots RN 1005	Quote Required
CAL81520	5 mL Model 1005 TLL Calibrated SYR NDL Sold Separately	Quote Required

Accessories

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.