

# 5 mL, Model 1005 SL SYR, Large Removable NDL, 22 ga, 2 in, point style 2

5 mL, Model 1005 SL SYR, Large Removable NDL, 22 ga, 2 in, point style 2



Part Number/REF: 81556

Unit Price: Quote Required

## Specifications

Volume	5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Sample Lock (SL)
Needle	Large RN
Gauge	22 gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	6.81 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.406 inches (10.3 mm)
Barrel Outer Diameter	0.530 inches (13.5 mm)
Sales Quantity	1 syringe per pack

## Accessories

Part Number/REF	Description	Unit Price
7779-01	22 gauge Large Hub RN NDL 2 in point style 2 6PK	Quote Required
30902	RN Nut	Quote Required
CAL81556	5 mL Model 1005 SL Calibrated SYR Large Removable NDL 22 ga 2 in point style 2	Quote Required

## Accessories



### Large Hub Removable Needles

These needles are available in 20-26 gauges. The standard length is 2 inches. Custom lengths are from 0.375 to 12 inches.

---

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.