

5 mL, Model 1005 C SYR, 1/4-28 Threads

5 mL, Model 1005 C SYR, 1/4-28 Threads



Part Number/REF: 81560

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	ChemSeal (C)
Needle	No Needle Available
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.406 inches (10.3 mm)
Barrel Outer Diameter	0.530 inches (13.5 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
13291	Plunger Assembly C 1005	Quote Required