

10 mL, Model 1010 LT SYR, NDL Sold Separately

10 mL, Model 1010 LT SYR, NDL Sold Separately



Part Number/REF: 81601

Unit Price: Quote Required

Specifications

| | |
|-------------------------|--------------------------|
| Volume | 10 mL |
| Syringe Type | Gastight |
| Syringe Series | 1000 |
| Termination | Luer Tip (LT) |
| Needle | Kel-F Hub |
| Gauge | No Needle Included |
| Point Style | No Needle Included |
| Needle Length | No Needle Included |
| Std. Needle Dead Volume | N/A |
| Fluid Path | Borosilicate glass, PTFE |
| Autoclavable | Yes to 115°C |
| Gas Sterilizable | Yes |
| Pressure | 200 psig (13.8 bar) |
| Maximum Temperature | 115°C / 239°F |
| Minimum Temperature | 10°C / 50°F |
| Barrel Inner Diameter | 0.574 inches (14.6 mm) |
| Barrel Outer Diameter | 0.695 inches (17.7 mm) |
| Sales Quantity | 1 syringe per pack |

Accessories

| Part Number/REF | Description | Unit Price |
|-----------------|--|----------------|
| 13231 | Plunger Assembly LTLTNTLL w wo Slots RN 1010 | Quote Required |
| CAL81601 | 10 mL Model 1010 LT Calibrated SYR NDL Sold Separately | Quote Required |

Accessories

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.