

10 mL, Model 1010 TLL Waters SYR, Waters HPLC Pump Priming Syringe

10 mL, Model 1010 TLL Waters SYR, Waters HPLC Pump Priming Syringe



Part Number/REF: 81610

Unit Price: Quote Required

Specifications

Volume	10 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Needle	Metal (N) or Kel-F
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.574 inches (14.6 mm)
Barrel Outer Diameter	0.695 inches (17.7 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
90049	25s gauge Metal Hub NDL 197 in point style 3 6PK	Quote Required
90149	25s gauge KelF Hub NDL 197 in point style 3 6PK	Quote Required

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.