

10 mL, Model 1010 LTSN SYR, Cemented NDL, Custom gauge, length, point style

10 mL, Model 1010 LTSN SYR, Cemented NDL, Custom gauge, length, point style



Part Number/REF: 81614
Unit Price: Quote Required

Custom Specifications

Custom Point Style	2
--------------------	---

Specifications

Custom Needle Length	0.75 inch (19.1 mm)
Volume	10 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip Cemented Needle Special (LTSN)
Needle	Cemented Needle
Std. Needle Dead Volume	N/A
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.574 inches (14.6 mm)
Barrel Outer Diameter	0.695 inches (17.7 mm)
Sales Quantity	1 syringe per pack
Custom Gauge	Custom Gauge
Custom Point Style	Custom Point Style

Accessories

Part Number/REF	Description	Unit Price
13231	Plunger Assembly LTLTNTLL w wo Slots RN 1010	Quote Required

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.