

10 mL, Model 1010 LTN SYR, Cemented NDL, 22 ga, 2 in, point style 3

10 mL, Model 1010 LTN SYR, Cemented NDL, 22 ga, 2 in, point style 3



Part Number/REF: 81616

Unit Price: Quote Required

Specifications

Volume	10 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip Cemented Needle (LTN)
Needle	Cemented Needle
Gauge	22 gauge
Point Style	3
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	6.81 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.574 inches (14.6 mm)
Barrel Outer Diameter	0.695 inches (17.7 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
13231	Plunger Assembly LTLTNTLL w wo Slots RN 1010	Quote Required

CAL81616

10 mL Model 1010 LTN Calibrated SYR Cemented NDL 22 ga 2 in point style 3

Quote Required

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.