

10 mL, Model 1010 RN SYR, Large Removable NDL, 22 ga, 2 in, point style 2

10 mL, Model 1010 RN SYR, Large Removable NDL, 22 ga, 2 in, point style 2



Part Number/REF: 81630
Unit Price: Quote Required

Specifications

Volume	10 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	22 gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0161 inches (0.41 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	6.81 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.574 inches (14.6 mm)
Barrel Outer Diameter	0.695 inches (17.7 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7779-01	22 gauge Large Hub RN NDL 2 in point style 2 6PK	Quote Required
13231	Plunger Assembly LTLTNTLL w wo Slots RN 1010	Quote Required
30902	RN Nut	Quote Required

CAL81630

10 mL Model 1010 RN Calibrated SYR Large Removable NDL 22 ga 2 in point style 2

Quote Required

Accessories



Large Hub Removable Needles

These needles are available in 20-26
gauges. The standard length is 2
inches. Custom lengths are from
0.375 to 12 inches.

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.