

25 μ L, Model 1802 N SYR, Cemented NDL, 22s ga, 2 in, point style 2

25 μ L, Model 1802 N SYR, Cemented NDL, 22s ga, 2 in, point style 2



Part Number/REF: 84878

Unit Price: Quote Required

Specifications

| | |
|-------------------------|-------------------------------|
| Volume | 25 μ L |
| Syringe Type | Gastight |
| Syringe Series | 1800 |
| Termination | Cemented Needle (N) |
| Needle | Cemented Needle |
| Gauge | 22s gauge |
| Point Style | 2 |
| Needle Length | 2 inch (51 mm) |
| Needle Inner Diameter | 0.0059 inches (0.15 mm) |
| Needle Outer Diameter | 0.0283 inches (0.72 mm) |
| Std. Needle Dead Volume | 1.13 μ L |
| Fluid Path | SST, borosilicate glass, PTFE |
| Autoclavable | No |
| Gas Sterilizable | Yes |
| Pressure | 1000 psig (68.9 bar) |
| Maximum Temperature | 50°C / 122°F |
| Minimum Temperature | 10°C / 50°F |
| Barrel Inner Diameter | 0.029 inches (0.729 mm) |
| Barrel Outer Diameter | 0.305 inches (7.75 mm) |
| Sales Quantity | 1 syringe per pack |

Accessories

| Part Number/REF | Description | Unit Price |
|-----------------|----------------------------|----------------|
| 1160-01 | Plunger Assembly N RN 1802 | Quote Required |
| 32135 | Holder Assembly | Quote Required |
| 32176 | Glass Barrel N 1802 | Quote Required |

CAL84878

25 microL I Model 1802 N Calibrated SYR Cemented NDL 22s ga 2 in point style 2

Quote Required

Accessories



Chaney Reproducibility Adapter

This adapter is built for reproducible injections. The adjustable stop prevents plunger damage and is easy to install in your laboratory.

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.