

100 µL, Model 1810 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2

100 µL, Model 1810 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2



Part Number/REF: 84886
Unit Price: Quote Required

Specifications

Volume	100 µL
Syringe Type	Gastight
Syringe Series	1800
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0059 inches (0.15 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	1.13 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.057 inches (1.46 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7758-03	22s gauge Small Hub RN NDL 2 in point style 2 6PK	Quote Required
1161-02	Plunger Assembly N RN 1810	Quote Required
30902	RN Nut	Quote Required

32135	Holder Assembly	Quote Required
32190	Glass Barrel RN 1810	Quote Required
CAL84886	100 microL a Model 1810 RN Calibrated SYR Small Removable NDL 22s ga 2 in point style 2	Quote Required

Accessories



Chaney Reproducibility Adapter

This adapter is built for reproducible injections. The adjustable stop prevents plunger damage and is easy to install in your laboratory.



Small Hub Removable Needles

These needles are available in 18-34 gauges. The standard length is 2 inches. Custom lengths are from 0.375 to 12 inches.

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.