

250 µL, Model 1825 RN SYR, Large Removable NDL, 22s ga, 2 in, point style 2

250 µL, Model 1825 RN SYR, Large Removable NDL, 22s ga, 2 in, point style 2



Part Number/REF: 84889

Unit Price: Quote Required

Specifications

Volume	250 µL
Syringe Type	Gastight
Syringe Series	1800
Termination	Removable Needle (RN)
Needle	Large RN
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.0059 inches (0.15 mm)
Needle Outer Diameter	0.0283 inches (0.72 mm)
Std. Needle Dead Volume	1.13 µL
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	500 psig (34.5 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.091 inches (2.30 mm)
Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
7779-03	22s gauge Large Hub RN NDL 2 in point style 2 6PK	Quote Required
30902	RN Nut	Quote Required
32135	Holder Assembly	Quote Required

32191	Glass Barrel RN 1825	Quote Required
4871-01	Plunger Assembly N 1825	Quote Required
CAL84889	250 microL y Model 1825 RN Calibrated SYR Large Removable NDL 22s ga 2 in point style 2	Quote Required

Accessories



Chaney Reproducibility Adapter

This adapter is built for reproducible injections. The adjustable stop prevents plunger damage and is easy to install in your laboratory.



Large Hub Removable Needles

These needles are available in 20-26 gauges. The standard length is 2 inches. Custom lengths are from 0.375 to 12 inches.

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.