

100 mL, Model 1100 TLL SYR, NDL Sold Separately

100 mL, Model 1100 TLL SYR, NDL Sold Separately



Part Number/REF: 86020

Unit Price: Quote Required

Specifications

Fluid Path	Borosilicate glass, PTFE, PCTFE
Volume	100 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	PTFE Luer Lock (TLL)
Needle	Metal (N) or Kel-F
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Autoclavable	No
Gas Sterilizable	Yes
Pressure	100 psig (6.9 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	1.28 inches (32.6 mm)
Barrel Outer Diameter	1.45 inches (36.9 mm)
Sales Quantity	1 syringe per pack

Accessories

Part Number/REF	Description	Unit Price
CAL86020	100 mL a Model 1100 TLL Calibrated SYR NDL Sold Separately	Quote Required

of Hamilton Company or its licensors in the U.S. and/or other countries.