

# 500 mL, Model S0500 Trachael SYR, No Needle Available

500 mL, Model S0500 Trachael SYR, No Needle Available



Part Number/REF: 86301

Unit Price: Quote Required

## Specifications

Volume	500 mL
Syringe Type	Gastight
Syringe Series	Super Syringe
Termination	Special
Needle	No Needle Available
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Buna-N, Nickel plated brass, PTFE, Silicone
Autoclavable	No
Gas Sterilizable	Yes
Pressure	30 psig (2.1 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	0°C / 32°F
Barrel Inner Diameter	2.44 inches (62.0 mm)
Barrel Outer Diameter	2.95 inches (74.9 mm)
Sales Quantity	1 syringe per pack

## Accessories

Part Number/REF	Description	Unit Price
18050	Super Syringe Barrel S0500	Quote Required

of Hamilton Company or its licensors in the U.S. and/or other countries.