

# 1 L, Model S1000 Trachael SYR, No Needle Available

---

1 L, Model S1000 Trachael SYR, No Needle Available



Part Number/REF: 86302

Unit Price: Quote Required

## Specifications

Volume	1 L
Syringe Type	Gastight
Syringe Series	Super Syringe
Termination	Special
Needle	No Needle Available
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Buna-N, Nickel plated brass, PTFE, Silicone
Autoclavable	No
Gas Sterilizable	Yes
Pressure	30 psig (2.1 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	0°C / 32°F
Barrel Inner Diameter	3.44 inches (87.4 mm)
Barrel Outer Diameter	3.95 inches (100 mm)
Sales Quantity	1 syringe per pack

---

© 2017 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.