

# 2 L, Model S2000 TLL SYR, NDL Sold Separately

2 L, Model S2000 TLL SYR, NDL Sold Separately



Part Number/REF: 86314

Unit Price: Quote Required

## Specifications

Volume	2 L
Syringe Type	Gastight
Syringe Series	Super Syringe
Termination	PTFE Luer Lock (TLL)
Needle	Metal (N) or Kel-F
Gauge	No Needle Included
Point Style	No Needle Included
Needle Length	No Needle Included
Std. Needle Dead Volume	N/A
Fluid Path	Buna-N, Nickel plated brass, PTFE, Silicone
Autoclavable	No
Gas Sterilizable	Yes
Pressure	30 psig (2.1 bar)
Maximum Temperature	80°C / 176°F
Minimum Temperature	0°C / 32°F
Barrel Inner Diameter	3.44 inches (87.4 mm)
Barrel Outer Diameter	3.95 inches (100 mm)
Sales Quantity	1 syringe per pack

## Accessories

Part Number/REF	Description	Unit Price
18026	TLL Connector S0500 S1000 S1500 S2000	Quote Required

