

2.5 mL, Model 1002 LTSN SYR, Cemented NDL, Custom gauge, length, point style

2.5 mL, Model 1002 LTSN SYR, Cemented NDL, Custom gauge, length, point style



Part Number/REF: 81414
Unit Price: Quote Required

Custom Specifications

Custom Point Style	2
--------------------	---

Specifications

Custom Needle Length	0.75 inch (19.1 mm)
Volume	2.5 mL
Syringe Type	Gastight
Syringe Series	1000
Termination	Luer Tip Cemented Needle Special (LTSN)
Needle	Cemented Needle
Std. Needle Dead Volume	N/A
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	200 psig (13.8 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.287 inches (7.29 mm)
Barrel Outer Diameter	0.407 inches (10.3 mm)
Sales Quantity	1 syringe per pack
Custom Gauge	Custom Gauge
Custom Point Style	Custom Point Style

Accessories

Part Number/REF	Description	Unit Price
1360-01	Plunger Assembly LTLTNTLL w amp wo Slots C 1002	Quote Required

Accessories



PB600-1 Repeating Dispenser

The PB600-1 is a sampling device designed for repeated dispensing of equal volumes. The instrument delivers 1/50th of the syringe volume.

© 2016 Hamilton Company. All rights reserved. Hamilton Company, product and/or service names referenced are either trademarks or registered trademarks of Hamilton Company or its licensors in the U.S. and/or other countries.