

Microliter 600 Series Syringes



The 600 series syringes are a half stroke version of the original 700 series syringes. This series of syringes is ideal for dispensing volumes from 0.25 µL up to 5 µL. The 600 series is often used in applications that require one handed operation because of the short stroke length and the reinforced plunger.

Downloads

Technical Information

Contact Us

+1-800-648-5950

sales@hamiltoncompany.com (<mailto:sales@hamiltoncompany.com>)

Displaying 1 - 4 of 4 products

Part Number/REF 	Description 	Unit Price 	Quantity
7632-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/25-microL-Model-62-RN-SYR-NDL-Sold-Separately)	2.5 µL, Model 62 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/25-microL-Model-62-RN-SYR-NDL-Sold-Separately)	Quote Required	1 Add to Cart +
87942 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/25-microL-Model-62-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-3)	2.5 µL, Model 62 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 3 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/25-microL-Model-62-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-3)	Quote Required	1 Add to Cart +
7633-01 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/5-microL-Model-65-RN-SYR-NDL-Sold-Separately)	5 µL, Model 65 RN SYR, NDL Sold Separately (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/5-microL-Model-65-RN-SYR-NDL-Sold-Separately)	Quote Required	1 Add to Cart +
87943 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/5-microL-Model-65-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-3)	5 µL, Model 65 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 3 (http://www.hamiltoncompany.com/products/syringes-and-needles/general-syringes/microliter-syringes/600-series/5-microL-Model-65-RN-SYR-Small-Removable-NDL-22s-ga-2-in-point-style-3)	Quote Required	1 Add to Cart +

Display

20 50 All